



Recombinant Escherichia coli Starvation-sensing protein rspA (rspA)

Product Code	CSB-EP327780ENV
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P38104
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli (strain K12)
Purity	>85% (SDS-PAGE)
Sequence	MKIVKAEV FV TCPGRNFVTL KITTEDGITG LGDATLNGRE LSVASYLQDH LCPQLIGRDA HRIEDIWQFF YKGAYWRRGP VTMSAISAVD MALWDIKAKA ANMPYQLLG GASREGVMVY CHTTGHSIDE ALDDYARHQE LGFKAIRVQC GIPGMKTTYG MSKGGKGLAYE PATKGQWPEE QLWSTEKYLD FMPKLFDAVR NKFGFNEHLL HDMHHLRTP I EAARFGKSIE DYRMFWMEDP TPAENQECFR LIRQHTVTPI AVGEVFNSIW DCKQLIEEQL IDYIRTTLTH AGGITGM RRI ADFASLYQVR TGSHGPS DLS PVCMAAALHF DLWVPNFGVQ EYMGYSEQML EVFPHNWTFD NGYMHPGDKP GLGIEFDEKL AAKYPYEPAY LPVARLEDGT LWNW
Source	E.coli
Target Names	rspA
Protein Names	Recommended name: Starvation-sensing protein rspA
Expression Region	1-404
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.