



Recombinant *Saccharomyces cerevisiae* Uncharacterized protein YBL086C (YBL086C)

Product Code	CSB-MP327792SVG
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P38177
Product Type	Recombinant Protein
Immunogen Species	<i>Saccharomyces cerevisiae</i> (strain ATCC 204508 / S288c) (Baker's yeast)
Purity	>85% (SDS-PAGE)
Sequence	MPFNHNSKAK RPKFLLDLQI KELVNIPQSS GYCYTKWRLK DGTGTSGHKV ALDGEHQTTT TQSRGTTKHV HVQHHRQAWN YSLDKPILVK LHLDKNGRFL KKILVLDVFF EFADANSSLT SSSSPNGKVK KTTYANATAL TATGNNSYSQ KITGKLLLGT VDIDITEYVK EDETPTTNRFL LKHSKVNSI INVSLQLKLV RGSYEDFNIS KSFTNGQLAN YRPGINTILD NTSELSSPTS TTNQMSPKNT FSNFNGIGTT VAKPGTNATG NSTSIKSPTS TNHKSSEMTT KPGLSTTISS SMSPLIESLY QKTFKLPWDP RPGEFTPREC VEDILQGGNG WAKNEKGINL IDLQALRLNE MEEYYNPNY GNNLGNKASS WPPNPSDDGY STMGKREYLE KKQNWSHMSR AQRALRTHN DEDNENTAND KGSDKDNNSV EDNNPTDFLT DRIREKNWS IITPSG
Source	Mammalian cell
Target Names	YBL086C
Protein Names	Recommended name: Uncharacterized protein YBL086C
Expression Region	1-466
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.