



Recombinant *Saccharomyces cerevisiae* Uncharacterized GTP-binding protein OLA1 (OLA1)

| | |
|--------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Code | CSB-EP327800SVG-B |
| Storage | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |
| Uniprot No. | P38219 |
| Product Type | Recombinant Protein |
| Immunogen Species | <i>Saccharomyces cerevisiae</i> (strain ATCC 204508 / S288c) (Baker's yeast) |
| Purity | >85% (SDS-PAGE) |
| Sequence | MPPKKQVEEK KVLLGRPGNN LKAGIVGLAN VGKSTFFQAI TRCPLGNPAN YPFATIDPEE ARVIVSPRF DKLCEIYKKT ASEVPAHLTV YDIAGLTKGA SAGEGLGNAF LSHIRSVDSI YQVVRCFDDA EIIHVEGDVD PVRDLEIINQ ELRLKDIEFA QKALEGAEKI AKRGGQSLEV KQKKEEMDLI TKIIKLLESG QRVANHSWTS KEVEIINSMF LLTAKPCIYL INLSERDYIR KKNKHLLRIK EWVDKYSPGD LIIPFSVSLE ERLSHMSPED AEEELKKLQT ISALPKIITT MRQKLDLISF FTCGPDEVRE WTIRRGTKAP QAAGVIHNDL MNTFILAQVM KCEDVFEYKD DSAIKAAGKL MQKGKDYVVE DGDIIYFRAG AGKN |
| Source | <i>E.coli</i> |
| Target Names | OLA1 |
| Protein Names | Recommended name: Uncharacterized GTP-binding protein OLA1 Alternative name(s): Obg-like ATPase 1 |
| Expression Region | 1-394 |
| Notes | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week. |
| Tag Info | Tag type will be determined during the manufacturing process. |
| Protein Length | full length protein |
| Reconstitution | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference. |
| Shelf Life | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |