



Recombinant Zea mays Dof zinc finger protein MNB1A (MNB1A)

Product Code	CSB-MP327857ZAX
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P38564
Product Type	Recombinant Protein
Immunogen Species	Zea mays (Maize)
Purity	>85% (SDS-PAGE)
Sequence	MQEASSAAAA GAEPGRRAAQ HQFAGVDLRR PKGYAAPAPA PAVGEGDPCP RCASRDTKFC YNNYNTSQP RHFCCKGRRY WTKGGTLRNV PVGGGTRKKP SSSSSSSSYV AAADADRQPK KKPASKKRRV VAPAPELATA ADPGKTATTT TTTSEITTET GALEDSDSLA HLLQPGTED AEVALGLGL SDFPSAGKAV LDEEDSFVWP AASFDMGACW AGAGFADPDP ACIFLNLP
Source	Mammalian cell
Target Names	MNB1A
Protein Names	Recommended name: Dof zinc finger protein MNB1A
Expression Region	1-238
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.