



Recombinant Escherichia coli Periplasmic protein torT (torT)

Product Code	CSB-BP327872ENV
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P38683
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli (strain K12)
Purity	≥85% (SDS-PAGE)
Sequence	AD NLLRWHAQH FTVQASTPLK AKRAWKLCAL YPSLKDSYWL SLNYGMQEAA RRYGVDLKVL EAGGYSQLAT QQAQIDQCKQ WGAEAILLGS STTSFPDLQK QVASLPVIEL VNAIDAPQVK SRVGVFWFQM GYQPGRYLVQ WAHGKPLNVL LMPGPDNAGG SKEMVEGFRA AIAGSPVRIV DIALGDNDIE IQRNLLQEML ERHPEIDVVA GTAIAAEAAM GEGRNKLTPL TVVSFYLSHQ VYRGLKRGRV IMAASDQMVW QGELAVEQAI RQLQGQSVSD NVSPILVLT PKNADREHIR RSLSPGGFRP VYFYQHTSAA KK
Source	Baculovirus
Target Names	torT
Protein Names	Recommended name: Periplasmic protein torT
Expression Region	19-342
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.