



Recombinant *Saccharomyces cerevisiae*

Uncharacterized protein YHR033W (YHR033W)

Product Code	CSB-YP327873SVG
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P38690
Product Type	Recombinant Protein
Immunogen Species	<i>Saccharomyces cerevisiae</i> (strain ATCC 204508 / S288c) (Baker's yeast)
Purity	≥85% (SDS-PAGE)
Sequence	MTKAYTIVIK LGSSSLVDES TKEPKLSIMT LIVETVTNLK RMGHKVIIVS SGGIAVGLDA LNIPHKPKQL SEVQAIAAVG QGRLIARWNM LFSQYGEQTA QILLTRNDIL RWNQYNNARN TINELLAMGV IPIVNENDTL SISEIEFGDN DTLSAITAAL VGADFLFLT DVDCLYTDNP RTNPDARPIV LPELSEGLP GVNTSSGSGS EVGTGGMRTK LIAADLASNA GIETIVMKSD RPEYVPKIVD YIQHFRPPR HIGNGTQQQF LDLQDTELEQ LRRYDVPMT KFLANDNKHK LKNREFWILH GLITKGAIID DENSYDKLLS KDMASLTPNA VIEVRDNFHE LECVDLKIGQ RLPNGELDIS KPIQSVGCVR SNYTSLELAK IKGLPSEKIH DVLGYSVSEY VAHRENIAFP PQF
Source	Yeast
Target Names	YHR033W
Protein Names	Recommended name: Uncharacterized protein YHR033W
Expression Region	1-423
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.