



Recombinant *Saccharomyces cerevisiae* Uncharacterized protein YHR003C (YHR003C)

Product Code	CSB-MP327885SVG
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P38756
Product Type	Recombinant Protein
Immunogen Species	<i>Saccharomyces cerevisiae</i> (strain ATCC 204508 / S288c) (Baker's yeast)
Purity	>85% (SDS-PAGE)
Sequence	MANNTWKLIA TTALISVFST QLAKSVWKEY KLSCAANKNK TVSRPRQYDD HLFREQLARN YAFLGEEGMR KIKEQYIVIV GAGEVGSWVC TMLIRSGCQK IMIIDPENIS IDSLNTHCCA VLSDIGKPKV QCLKEHLSKI APWSEIKARA KAWTKENSHD LIFADGESPT FIVDCLDNLE SKVDLLEYAH HNKIDVISSM GVATKSDPTR VSINDISMTE FDPISRCVRR KLRKRGATG ISVVFSSNEML DPRRDDILSP IDCEHRAINA VRDEALRHLP ELGTMPGIFG LSIATWILTK VSGYPMKENE VKNRLKFYDS ILETQKQMA RLNENKERSS LLGLEEVGYI VEEMFRGKSP ISGYSTKLAL TKWEANKEIS LTNVVLMTKE EQEIHEKRIL LDGEKLTAVY SEEVLDIFIER LFKEEEYYSS
Source	Mammalian cell
Target Names	TCD1
Protein Names	Recommended name: Uncharacterized protein YHR003C
Expression Region	1-429
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.