



Recombinant *Saccharomyces cerevisiae* Pre-rRNA-processing protein IPI1 (IPI1)

Product Code	CSB-MP327892SVG
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P38803
Product Type	Recombinant Protein
Immunogen Species	<i>Saccharomyces cerevisiae</i> (strain ATCC 204508 / S288c) (Baker's yeast)
Purity	≥85% (SDS-PAGE)
Sequence	MTKSRKQKQK KQDFLRKKLK VGKPKEKARN ATDTSFVSKT ISIRNQHLDQ NPHDLTKRLT LLKHHNINVR KETLTTFQKS IPSIIKSRLM TPLLTQSIPL ICDESQQVRQ GLIDLVDIEIG SHDAEILKLH CNIFVLYINM AMTHIVTQIQ ADSTKFLSHL LKYCGDEVVR KSWVKLLNGV FGVLGWGQVG KNDSASIVQT KKRNAKYVTI HLNALYTLVE YGCQDERARS DGDTAETTED SGTLRNPYLI PDYPQPFEHL KLFTRCLKVQ DATSSGVNAT LLSLATQDID TRKAVFIEQF LPIVRKKIEV IIEGEGCGK SANKLKTLA KIFD
Source	Mammalian cell
Target Names	IPI1
Protein Names	Recommended name: Pre-rRNA-processing protein IPI1 Alternative name(s): Involved in processing IST2 protein 1
Expression Region	1-334
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.