



Recombinant *Saccharomyces cerevisiae* Chromatin modification-related protein YNG2 (YNG2)

Product Code	CSB-EP327893SVG
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P38806
Product Type	Recombinant Protein
Immunogen Species	<i>Saccharomyces cerevisiae</i> (strain ATCC 204508 / S288c) (Baker's yeast)
Purity	>85% (SDS-PAGE)
Sequence	MDPSLVLEQT IQDVSNLPSE FRYLLEEIGS NDCLKLIEEKK KYEQKESQIH KFIRQQGSIP KHPQEDGLDK EIKESLLKCQ SLQREKCVLA NTALFLIARH LNKLEKNIAL LEEDGVLAPV EEDGDMDSAA EASRESSVVS NSSVKKRRAA SSSGSVPPTL KKKKTSRTSK LQNEIDVSSR EKSVTPVSPS IEKKIARTKE FKNSRNGKGQ NGSPENEEED KTLYCFCQRV SFGEMVACDG PNCKYEFWHY DCVNLKEPPK GTWYCPECKI EMEKNLKRK RN
Source	E.coli
Target Names	YNG2
Protein Names	Recommended name: Chromatin modification-related protein YNG2 Alternative name(s): ESA1-associated factor 4 ING1 homolog 2
Expression Region	1-282
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.