



Recombinant *Saccharomyces cerevisiae* Uncharacterized isomerase YHR210C (YHR210C)

Product Code	CSB-MP327908SVG
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P38893
Product Type	Recombinant Protein
Immunogen Species	<i>Saccharomyces cerevisiae</i> (strain ATCC 204508 / S288c) (Baker's yeast)
Purity	>85% (SDS-PAGE)
Sequence	MSNNKAGGEY EVITIGDAKK LQATISELGA TLLDLKVNNE SIVLGYPDIH GYISDGYNYI GATVGRYANR IYKGMFSMED GPHQLTVNNC GNTNHSSISS FHLKKYKASK VQNPLDDLYI VEFTLLDDRT LPNEFPGDLA VNLKYTLNVA DMTLDLEYEA KLVSGEATPI NMTNHTYFNL NKTMDKKSIS GTEVRLCSDK SLEVSEGALI PTGKIVQRKI ATFDSSKPTI LQDDGPIYDY AFIVDENKNL KTTDSVSVNK LVPAFKAYHP ASRLSLEVST TEPTVLFYTG DNLCDGFTPR SGFAVEQGRY VDAINRDGWR DCVLLRRGEV YTSKTRYRFA V
Source	Mammalian cell
Target Names	YHR210C
Protein Names	Recommended name: Uncharacterized isomerase YHR210C EC= 5.-.-.-
Expression Region	1-341
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.