



# Recombinant *Saccharomyces cerevisiae* Protein CAJ1 (CAJ1)

|                          |                                                                                                                                                                                                                                                                                                                                                                                                                                                                  |
|--------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Product Code</b>      | CSB-YP327933SVG                                                                                                                                                                                                                                                                                                                                                                                                                                                  |
| <b>Storage</b>           | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.<br>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.                                                                                                                                                                               |
| <b>Uniprot No.</b>       | P39101                                                                                                                                                                                                                                                                                                                                                                                                                                                           |
| <b>Product Type</b>      | Recombinant Protein                                                                                                                                                                                                                                                                                                                                                                                                                                              |
| <b>Immunogen Species</b> | <i>Saccharomyces cerevisiae</i> (strain ATCC 204508 / S288c) (Baker's yeast)                                                                                                                                                                                                                                                                                                                                                                                     |
| <b>Purity</b>            | >85% (SDS-PAGE)                                                                                                                                                                                                                                                                                                                                                                                                                                                  |
| <b>Sequence</b>          | MVKETEYYDI LGIKPEATPT EIKKAYRRKA METHPKHPD DPDAQAKFQA<br>VGEAYQVLSD PGLRSKYDQF GKEDAVPQQG FEDASEYFTA IFGGDGFKDW<br>IGEFSLFKEL NEATEMFGKE DEEGTAATET EKADESTDGG MVKHDTNKA<br>SLKKDKLSKE QREKLMEMEK KRREDMMKQV DELAEKLNEK ISRYLIAVKS<br>NNLEEFTRKL DQEIEDLKL SFGLELLYLL ARVYKTKANN FIMSKKTYGI<br>SKIFTGTRDN ARSVKSAYNL LSTGLEAQA MEKMSEVNTD ELDQYERAKF<br>ESTMAGKALG VMWAMSKFEL ERKLKDV CNK ILNDKKVPSK ERIAKAKAML<br>FIAHKFASAR RSPEEAEEAR VFEELILGEQ EKEHKKHTVA R |
| <b>Source</b>            | Yeast                                                                                                                                                                                                                                                                                                                                                                                                                                                            |
| <b>Target Names</b>      | CAJ1                                                                                                                                                                                                                                                                                                                                                                                                                                                             |
| <b>Protein Names</b>     | Recommended name: Protein CAJ1                                                                                                                                                                                                                                                                                                                                                                                                                                   |
| <b>Expression Region</b> | 1-391                                                                                                                                                                                                                                                                                                                                                                                                                                                            |
| <b>Notes</b>             | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.                                                                                                                                                                                                                                                                                                                                                              |
| <b>Tag Info</b>          | Tag type will be determined during the manufacturing process.                                                                                                                                                                                                                                                                                                                                                                                                    |
| <b>Protein Length</b>    | full length protein                                                                                                                                                                                                                                                                                                                                                                                                                                              |
| <b>Reconstitution</b>    | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.                                                      |
| <b>Shelf Life</b>        | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.<br>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.                                                                                                                                                                               |