



Recombinant Escherichia coli D-mannonate oxidoreductase (uxuB)

Product Code	CSB-EP327942ENV-B
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P39160
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli (strain K12)
Purity	>85% (SDS-PAGE)
Sequence	MTTIVDSNLP VARPSWDHSR LESRIVHLGC GAFHRAHQAL YTHHLLSTD SDWGICEVNL MPGNDRVLIE NLKKQQLLYT VAEKGAESTE LKIIGSMKEA LHPEIDGCEG ILNAMARPQT AIVSLTVTEK GYCADAASGQ LDLNNPLIKH DLENPTAPKS AIGYIVEALR LRREKGLKAF TVMSCDNVRE NGHVAKVAVL GLAQARDPQL AAWIEENVTF PCTMVDRIVP AATPETLQEI ADQLGVYDPC AIACEPFRQW VIEDNFVNGR PDWDKVGAFQ VADVVPFEMM KLRMLNGSHS FLAYLGYLGG YETIADTVTN PAYRKAALFAL MMQEQAPTLS MPEGTDLNAY ATLLIERFSN PSLRHRTWQI AMDGSQKLPQ RLLDPVRLHL QNGGSRHLA LGVAGWMRYT QGVDEQGNAI DVVDPMLAEF QKINAQYQGA DRVKALLGLS GIFADDLPQN ADFVGAVTAA YQQLCERGAR ECVAAL
Source	E.coli
Target Names	uxuB
Protein Names	Recommended name: D-mannonate oxidoreductase EC= 1.1.1.57 Alternative name(s): Fructuronate reductase
Expression Region	1-486
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.