



Recombinant Human UPF0533 protein C5orf44 (C5orf44)

Product Code	CSB-BP004022HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	A5PLN9
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MEVNPPKQEH LLALKVMRLT KPTLFTNIPV TCEEKDLPGD LFNQLMRDDP STVNGAEVLM LGEMLTLPQN FGNIFLGETF SSYISVHNSD NQVVKDILVK ADLQTSSQRL NLSASNAAVA ELKPDCCIDD VIHHEVKEIG THILVCAVSY TTQAGEKMYF RKFFKFQVLK PLDVKTKFYN AESDLSSVTD EVFLEAQIQN MTTSPMFMEK VSLEPSIMYN VTELNSVSQA GECVSTFGSR AYLQPM DTRQ YLYCLKPKNE FAEKAGIIG VTVIGKLDIV WKTNLGERGR LQTSQLQRMA PGYGDVRLSL EAIPDTVNLE EPFHITCKIT NCSERTMDLV LEMCNTNSIH WCGISGRQLG KLHPSSSLCL ALTLLSSVQG LQSISGLRLT DTFLKRTY EY DDIAQVCVVS SAIKVES
Source	Baculovirus
Target Names	TRAPPC13
Protein Names	Recommended name: UPF0533 protein C5orf44
Expression Region	1-417
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.