



Recombinant Human Protein FAM149A (FAM149A)

Product Code	CSB-EP008082HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	A5PLN7
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	<p>MKAAVLDLGS LLAKLFETST APPAGPSSRP SGGAAAAGSG GSRAGTPLGT APTLLRALAP DSPSASRRSP APLLSSPYSR GSAASRAAGA VGTLISWPSS PRAGKAPPQP PTPSGGGCSP ARLVVPARPP SGPGGVWAAL PRNPLQPGPG ERELGACVAP GAGPRTLFLT LPDIGEEGAS DGDSGDGEAR GLSEGRRRRHG FTVRSKDSL P THFTRNVQKA IDKYTCKSL S SFSSSGSHTP TGAHTSWSGS ATQSSTTGSS TERGSVYSWR DDEFDEASSQ SVQRLLWEVE EMLFEGKVN P QTQSLLAECG EWTRRSLHLR VLGRQLILPT DKGVQHFQGS TPASAVHRPP LSACGHSSNI RELCISGSQI VPAALSASAL PGPDDTGVAD LTARSSLEEE VYHVDGKIEE YFAFDRKEDD DECLEQKPAQ PGRKWRKLGL PPVSPRDCVK DAVAAEVFDH VWTNMVELLE ELIRKHWETT LTEGKKQRET LKVAGNRFP H VLVPHAHADG ASGPPSGHAE AHGISLASRL NPPQIHFFSS SFYSDMNGVM TIQAKPLQRR PAYFADRTQN EKEDKASGGG AGALSSAPHR LGRASDTHGL SPSAKKTPVP WRLPSLASDS QRLKTPNIYS DEVLRGTKLP TGVDHMASPL VQTSRSRFP P LVTETRGQNT AVPGCRLVSY RGRHLQNRVL SAMPDGT ERS RLRERTATLE RLSRPSTTHT FRQSDTPRKS SLTQMEFAAH TWTGQSILTG SQYVPKSFQR TTLTFKRRFQ VTS</p>
Source	E.coli
Target Names	FAM149A
Protein Names	Recommended name: Protein FAM149A
Expression Region	1-773
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.