



Recombinant *Synechococcus elongatus* Bicarbonate-binding protein CmpA (cmpA)

Product Code	CSB-EP328010FPY
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P39660
Product Type	Recombinant Protein
Immunogen Species	<i>Synechococcus elongatus</i> (strain PCC 7942) (<i>Anacystis nidulans</i> R2)
Purity	>85% (SDS-PAGE)
Sequence	GGGT SSTTQPTAAG ASDLEVKTIK LGYIPIFEAA PLIIGREKGF FAKYGLDVEV SKQASWAAAR DNVILGSAGG GIDGGQWQMP MPALLTEGAI SNGQKVPMYV LACLSTQGNG IAVSNQLKAQ NLGLKLAPNR DFILNYPQTS GRKFKASYTF PNANQDFWIR YWFAAGGIDP DKDIELLTVP SAETLQNMNRN GTIDCFSTGD PWPSRIAKDD IGYQAALTGQ MWPYHPEEFL ALRADWVDKH PKATLALLMG LMEAQQWCDQ KANRAEMAKI LSGRNFFNVP VSILQPILEG QIKVGADGKD LNNFDAGPLF WKSPRGSVSY PYKGLTLWFL VESIRWGFNK QVLPDIAAAQ KLNDRVTRED LWQEAAKKG VPAADIPTGS TRGTETFFDG ITYNPDSPQA YLQSLKIKRA
Source	E.coli
Target Names	cmpA
Protein Names	Recommended name: Bicarbonate-binding protein CmpA
Expression Region	37-450
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.