



Recombinant *Synechococcus elongatus* Protein sphX (sphX)

Product Code	CSB-EP328011FPY-B
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P39665
Product Type	Recombinant Protein
Immunogen Species	<i>Synechococcus elongatus</i> (strain PCC 7942) (<i>Anacystis nidulans</i> R2)
Purity	≥85% (SDS-PAGE)
Sequence	SAQRALVTAD GSSTVFPISE AVAAEFQKRN KNINVTVGVS GTGGGFKRFC NGEIDIANAS RPIKKEEVEA CRKKGIRYIE LPVAFDALTV VVNKSNPVNS ITTAELAKIF GRDAEKTTN WRQVKSSFPN LPLRVYAPGT DSGTYDYFNE AILNKKGTRG DLTASEDDNI LVQGVSRDRG GIGFFGFSYY EENKGKALKAL AVVNSNGKAV MPSVQNVLNG TYDPLARPVF IYVSEQAACK ANVRSFVNFY LQNAGKLSRE VGFVPLPAKA YTAATQRFRS NKTGTVFAGK SLVGGSIEDL LKAEGIN
Source	E.coli
Target Names	sphX
Protein Names	Recommended name: Protein sphX
Expression Region	31-337
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.