



Recombinant *Saccharomyces cerevisiae* Putative metallophosphoesterase YNL217W (YNL217W)

Product Code	CSB-MP328080SVG
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P40152
Product Type	Recombinant Protein
Immunogen Species	<i>Saccharomyces cerevisiae</i> (strain ATCC 204508 / S288c) (Baker's yeast)
Purity	>85% (SDS-PAGE)
Sequence	LYIF KFDNPRLSPP VSLLP TISTL KKIEHVTDLN KEYV FVG DVH GNYDEFIELI DDKIGGLGEN ITMILLGDFI HKGPDS DKVV SYILNHK DQV KCVLGNHEIL VMMAYLNP DF SKWVRRPKLM TPLTFSTETN FIPQDISKIS NAHGRLAREL GFSKLSQLAE HCSMAIELDL DITGDILFGA HAGMVP GDFM KPNQIPGVSS LSNMKYVDKK NWSKTSREKE NKNYVRWYTL WDKYGDHFSN AKVFGHDAS MGLNLRRQTK GLDTACIKNN LLSSMKVKYD IKKGQYDYEL IQVQCS
Source	Mammalian cell
Target Names	YNL217W
Protein Names	Recommended name: Putative metallophosphoesterase YNL217W EC= 3.1.-.-
Expression Region	27-326
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.