



Recombinant *Saccharomyces cerevisiae* Anaphase-promoting complex subunit MND2 (MND2)

Product Code	CSB-YP328139SVG
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P40577
Product Type	Recombinant Protein
Immunogen Species	<i>Saccharomyces cerevisiae</i> (strain ATCC 204508 / S288c) (Baker's yeast)
Purity	>85% (SDS-PAGE)
Sequence	MARALRDISL FNDIRKQNS AGAKHERYNM RDLRSKKNQH VNGIDDYEDD SLDRFIRRK SRVVKYIPSL SAYNVFNEFP YYPTSASQLL DGKLDEFMLL SEQYKSRLPK IRKLGWNRFK PIGINKTMYE LEMLSRARA QNAEGNNEED FRQHDSREED PRNNGSIGRV ILPHILQENE EYDTGEGVTG LHSM PND SMA ILANNSANNS QNEEVSEEDE ISYDYDAEFD HVVDEDDNEE GEVPGEGVEG IEVQRERIVP DDLLMRPTSL SRSLQQFVEE AHHLDRNPYD IDSDNDGEDS KVELDMNPFD EDDVGREHDY NSEYSQEPTS YGGITPD LAS NWRNWTRERI TSLDELMERR ARQQRGQD
Source	Yeast
Target Names	MND2
Protein Names	Recommended name: Anaphase-promoting complex subunit MND2 Alternative name(s): Meiotic nuclear division protein 2
Expression Region	1-368
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.