



# Recombinant *Drosophila melanogaster* Ornithine decarboxylase 2 (Odc2)

<b>Product Code</b>	CSB-MP328172DLU
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	P40808
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	<i>Drosophila melanogaster</i> (Fruit fly)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MVNGDLRIQY YDEELNIRKV IEQADLEHLD QALNICDLSS LERKLRLWHK LMPRIEPHYA VKCNDDPVVV KFLADLGTGF DCASKNELKL VLGLGVSPER IIFAHPCRPA SHLRYAKEQQ VVNGTVDNEY EIYKLRKHYP DSNLIVRFKS EAKKALCPLG DKYGCDAEAD AAALMLLAKA LGLKVTGTSF HVGSGCSEVE AYDRAIEKAE NIFKVGEMIG HKMELLDVGG GFPIDDEMFE EIAQSVNTS VELRFPDKRI RIIEPGRFF VEAAYTLICK VHAKREVRSK DGKLDTMMYY LNDGIFGAGA GMFYYPEEVA PELYLDEAES LPKLKSVIWG PSCDAMDKIS DLLLPNLNPG DLLGFRNMGA YTMPIASPFN GFDVPETRFF KAK
<b>Source</b>	Mammalian cell
<b>Target Names</b>	Odc2
<b>Protein Names</b>	Recommended name: Ornithine decarboxylase 2 Short name= ODC EC= 4.1.1.17 Alternative name(s): dODC2
<b>Expression Region</b>	1-393
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.