



Recombinant Candida albicans 40S ribosomal protein S1 (RPS1)

Product Code	CSB-EP328191CZD-B
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P40910
Product Type	Recombinant Protein
Immunogen Species	Candida albicans (strain SC5314 / ATCC MYA-2876) (Yeast)
Purity	>85% (SDS-PAGE)
Sequence	AVGKNKRLS KGKKGLKKKV VDPFTRKDFW DIKAPTTFEN RNVGKTLINR STGLKNAADG LKGRVFEVCL ADLQGSSEHS YRKIKLRVDE VQGKLLTNF HGLDFTSDKL RSLVRKWQSL VEANVTVKTS DDYVLRVFAI AFTKRQPNQI KKTYYAQSSK LREVRKKMIE IMQREVSNT LAQLTSKLIP EVIGREIEKS TQTIFPLQNV HIRKVKLLKQ PKFDLGLSLA LHGEGSTEEK GKKVSSGFKD VVLESV
Source	E.coli
Target Names	RPS1
Protein Names	Recommended name: 40S ribosomal protein S1 Alternative name(s): S3aE
Expression Region	2-256
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.