



Recombinant *Saccharomyces cerevisiae* Ribonuclease MRP protein subunit SNM1 (SNM1)

Product Code	CSB-MP328201SVG
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P40993
Product Type	Recombinant Protein
Immunogen Species	<i>Saccharomyces cerevisiae</i> (strain ATCC 204508 / S288c) (Baker's yeast)
Purity	>85% (SDS-PAGE)
Sequence	MNKDQAEKYQ ERS LRQKYNL LHVLP TLNSR ALSGLYYKNF HNSVKRYQIM LPEQLKSGKF CSHCGCVYVP NFNASLQLTT NTEQGDSDEL GGESMEGPKK CIQVNCLNCE KSKLFEWKSE FVVPTFGQDV SPMINSTSSG KVSYAVKKPQ KSKTSTGKER SKKRKLNSLT NLLSKRNQEK KMEKKKSSSL SLESFMKS
Source	Mammalian cell
Target Names	SNM1
Protein Names	Recommended name: Ribonuclease MRP protein subunit SNM1 Alternative name(s): RNA-processing protein SNM1 RNase MRP 22.5 kDa subunit
Expression Region	1-198
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.