



Recombinant *Drosophila melanogaster* Glutathione S-transferase S1 (GstS1)

| | |
|--------------------------|---|
| Product Code | CSB-EP328209DLU-B |
| Storage | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |
| Uniprot No. | P41043 |
| Product Type | Recombinant Protein |
| Immunogen Species | <i>Drosophila melanogaster</i> (Fruit fly) |
| Purity | >85% (SDS-PAGE) |
| Sequence | MADEAQAPPA EGAPPAEGEA PPPAEGAEGA VEGGEAAPPA EPAEPIKHSY TLFYFNVKAL AEPLRYLFAY GNQEYEDVRV TRDEWPALKP TMPMGQMPVL EVDGKRVHQS ISMARFLAKT VGLCGATPWE DLQIDIVVDT INDFRLKIAV VSYEPEDAIK EKKLVTLNAE VIPFYLEKLE QTVKDNDGHL ALGKLTWADV YFAGITDYMN YMVKRDLEP YPALRGVVDA VNALEPIKAW IEKRPVTEV |
| Source | <i>E. coli</i> |
| Target Names | GstS1 |
| Protein Names | Recommended name: Glutathione S-transferase S1 EC= 2.5.1.18 Alternative name(s): GST class-sigma 1 Glutathione S-transferase 2 |
| Expression Region | 1-249 |
| Notes | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week. |
| Tag Info | Tag type will be determined during the manufacturing process. |
| Protein Length | full length protein |
| Reconstitution | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference. |
| Shelf Life | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |