



Recombinant Transforming protein cbl (V-CBL)

Product Code	CSB-BP328377CEJ
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P23092
Product Type	Recombinant Protein
Immunogen Species	Cas-NS-1 murine leukemia virus
Purity	>85% (SDS-PAGE)
Sequence	ASAGGGCRRG PSFSPGSIPS LAAERAPDPP LAMAGNVKKS SGAGGGGSGG SGAGGLIGLM KDAFQPHHHH HHLSPHPPCT VDKMMVEKCW KLMDKVVRRLC QNPVALKNS PPYILDLLPD TYQHLRTVLS RYEGKMETLG ENEYFRVFME NLMKKTQTI SLFKEGKERM YEENSQPRRN LTKLSLIFSH MLAELKGIFP SGLFQGDTFR ITKADAAEFW RKAFGEKTIV PWKSFRQALH EVHPISSGLE AMALKSTIDL TCNDYISVFE FDIFTRLFQP WSSLLRNWNS LAVTHPGYMA FLTYDEVKAR LQKFIHKPGS YIFRLSCTRL GQWAIGYVTA DGNILQTIPH NKPLFQALID GFREGFYLFP DGRNQNPDLT GLCEPTPHFS
Source	Baculovirus
Target Names	V-CBL
Protein Names	Recommended name: Transforming protein cbl
Expression Region	1-390
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.