



Recombinant Aminopeptidase T

Product Code	CSB-MP328396TNN
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P23341
Product Type	Recombinant Protein
Immunogen Species	Thermus aquaticus
Purity	>85% (SDS-PAGE)
Sequence	MDAFTENLNK LAELAIRVGL NLEEGQEIVA TAPIEAVDFV RLLAEKAYEN GASLFTVLYG DNLIARKRLA LVPEAHLDRAPAWLYEGMAK AFHEGAARLA VSGNDPKALE GLPPERVGRA QQAQSRAYRP TLAITEFVT NWTIVPFAHP GWAKAVFPGL PEEEAQRLW QAIFQATRVD QEDPVAAWEA HNRVLHAKVA FLNEKRFHAL HFQGPGLDGLVGLAEGHLWQ GGATPTKKGR LCNPNLPTEE VFTAPHRERV EGVVRASRPL ALSGQLVEGL WARFEGGVAV EVGAEKGEEV LKKLLDTDEG ARRLGEVALV PADNPIAKTG LVFFDTLFDENAAASHIAFGQ AYAENLEGRP SGEEFRRRGG NESMVHVDWM IGSEEVDVDG LLEDGTRVPL MRRGRWVI
Source	Mammalian cell
Protein Names	Recommended name: Aminopeptidase T Short name= AP-T EC= 3.4.11.- Alternative name(s): Heat-stable aminopeptidase
Expression Region	1-408
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.