



Recombinant Agrobacterium tumefaciens Chloramphenicol acetyltransferase (cat)

Product Code	CSB-MP328398AYS
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P23364
Product Type	Recombinant Protein
Immunogen Species	Agrobacterium fabrum (strain C58 / ATCC 33970) (Agrobacterium tumefaciens (strain C58))
Purity	>85% (SDS-PAGE)
Sequence	MENYFESPFR GITLDKQVKS PNLVVGKYSY YSGYYHGHFS EDCARYLLPD EGADRLVIGS FCSIGSGAAF IMAGNQGHRN EWISTFPFFF MPEVPEFENA ANGYLPAGDT VIGNDVWIGS EAIIMPGITV GDGAVIGTRA LVTKDVEPYA IVGGNPAKTI RKRFD DDSIA LLEMKWWGW PAERLKAAMP LMTSGNVAAL YRFWRSDSL
Source	Mammalian cell
Target Names	cat
Protein Names	Recommended name: Chloramphenicol acetyltransferase EC= 2.3.1.28
Expression Region	1-209
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.