



Recombinant Schizosaccharomyces pombe Minor serine/threonine-protein phosphatase PP2A-1 catalytic subunit (ppa1)

Product Code	CSB-YP328418SXV
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P23635
Product Type	Recombinant Protein
Immunogen Species	Schizosaccharomyces pombe (strain 972 / ATCC 24843) (Fission yeast)
Purity	>85% (SDS-PAGE)
Sequence	MSVSGKIGEV DRWIEQLSRC EPLSEEDVIQ MCDLAKEVLS VESNVQSVRC PVTVCGDIHG QFHDLMELFN IGGPSPDTNY LFMGDYVDRG YHSVETVSL IAFKIRYPQR ITILRGNHES RQITQVYGFY DECLRKYGNA NVWQYFTDLF DYLPLTALIE DRIFCLHGGL SPSIDTLDHV RILDRVQVEP HEGPICDLLW SDPDDRPGWG ISPRGAGYTF GPDIAEAFNH NNGLDLIARA HQLVMEGYNW TTNHNVTIF SAPNYCYRCG NQAAIMGIDD HINYAFIQYD TAPRKEELHV TRRTPDYFL
Source	Yeast
Target Names	ppa1
Protein Names	Recommended name: Minor serine/threonine-protein phosphatase PP2A-1 catalytic subunit EC= 3.1.3.16
Expression Region	1-309
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.