



Recombinant Human Phosphoglycolate phosphatase (PGP)

Product Code	CSB-EP017873HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	A6NDG6
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MAAAEAGGDD ARCVRLSAER AQALLADVDT LLFDCDGV LW RGETAVPGAP EALRALRARG KRLGFITNNS SKTRAAYA EK LRRLGFGGPA GPGASLEVFG TAYCTALYLR QRLAGAPAPK AYVLGSPALA AELEAVGVAS VGVGPEPLQG EGPGDWLHAP LEPDVRAVVV GFDPHFSYMK LTKALRYLQQ PGCLLVGTNM DNRLPLENGR FIAGTGCLVR AVEMAAQRQA DIIGKPSRFI FDCVSQEYGI NPERTVMVGD RLDTDILLGA TCGLKTILT L TGVSTLGDVK NNQESDCVSK KKMVPDFYVD SIADLLPALQ G
Source	E.coli
Target Names	PGP
Protein Names	Recommended name: Phosphoglycolate phosphatase Short name= PGP Short name= PGPase EC= 3.1.3.18
Expression Region	1-321
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.