




Recombinant Vaccinia virus Protein OPG190 (OPG190), partial

Product Code	CSB-YP328444VAG
Abbreviation	Recombinant Vaccinia virus PS/HR protein, partial
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P24083
Storage Buffer	Tris-based buffer,50% glycerol
Product Type	Recombinant Proteins
Immunogen Species	Vaccinia virus (strain Lister) (VACV)
Purity	≥ 90% as determined by SDS-PAGE.
Sequence	YSTCTVPTMNNAKLTSTETSFNDKQKVTFCTDQGYHSSDPNAV CETDKWKYE NPCKKMCTVSDYVSELYDKPLYEVNSTMTLSCNGETKYFRCEEKNGNTSWN DTVTCPNAECQPLQLEHGSCQPVKEKYSFGHEYITINCDVG YEYIGASYISCTAN SWNVIPSCQQKCDMPSLSNGLISGSTFSIGGVIHL SCKSGFTLTGSPSSTCIDG KWNPIPTCVRSNEKFDPVDDGPDDDETDL SKLSKDVVQYEQEIESLEATYH
Research Area	Others
Source	Yeast
Target Names	OPG190
Protein Names	Recommended name: Plaque-size/host range protein Alternative name(s): Protein B5
Expression Region	18-279aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 6xHis-tagged
Mol. Weight	31.1 kDa
Protein Length	Partial
Image	



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.