



Recombinant Danio rerio Protein Wnt-1 (wnt1)

Product Code	CSB-YP328457DIL
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P24257
Product Type	Recombinant Protein
Immunogen Species	Danio rerio (Zebrafish) (Brachydanio rerio)
Purity	≥85% (SDS-PAGE)
Sequence	S LTGTGAVNNS GRWWGIVNVA SSGNLLTNSK NVQLVLDPSL ALLSRRQRKL IRQNPGILHA IAAGLHTAIK ECKWQFRNRR WNCPTTHSPN VFGKIVNRGC RETAFVFAIT SAGVTHAVAR SCSEGAIESC TCDYRRRGPG GPDWHWGGCS DNVEFGRMFG REFVDSSERG RDLRYLTNLH NNEAGRMTVA SEMQQECKCH GMSGSCTVRT CWMRLPSFRL VGDYLKDRFD GASRVVYANK GSNRASHRAD PRHLEPENPA HKLPSSRDV YFEKSPNFCS YNGKTGTHGT SGRTCNSSSP ALDGCELLCC GRGYKTRMEQ VTERCHCTFH WCCHVSLCNC TSTQTVHQCL
Source	Yeast
Target Names	wnt1
Protein Names	Recommended name: Protein Wnt-1
Expression Region	20-370
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.