



# Recombinant Type-2 restriction enzyme Accl (acclR)

<b>Product Code</b>	CSB-EP328478AWR-B
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	P24546
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Acinetobacter calcoaceticus
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	MDYYDRIREL TKNVPVELVD FEQPRDLART PTQASSNFIT NKEQGDWAED LVTRAINENS KNFVAVKYGK SDNLVAGENG FDTFYQDFQT ELDTIGKRPD LLIFKKTDFD TTLGFDVSI PHHQITDYVK KAIAGIEVRS SAFLIDKYEE AMQVRTQRFT EIAFQTRDKI LAEFLDVLDP PSRSKYITLL NTLTLETISI FDFKVPGWRS NERLIEVNNL FKRLKVAIKE IQKRDYLSIT PKVEDIKVVY KIETFNVPY FYFQVFFDKV YGISFEQILT IISNSDNDGV IFSVEKDQVN QNKTTIKINS KTGYPISKV DEPTHEIRK EMDRGRLLFY VTFKGGTAYL DLNLRITLG IEEAEF
<b>Source</b>	E.coli
<b>Target Names</b>	acclR
<b>Protein Names</b>	Recommended name: Type-2 restriction enzyme Accl Short name= R.Accl EC= 3.1.21.4 Alternative name(s): Endonuclease Accl Type II restriction enzyme Accl
<b>Expression Region</b>	1-366
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.