



Recombinant *Saccharomyces cerevisiae* PHO85 cyclin-1 (PCL1)

Product Code	CSB-MP328505SVG
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P24867
Product Type	Recombinant Protein
Immunogen Species	<i>Saccharomyces cerevisiae</i> (strain ATCC 204508 / S288c) (Baker's yeast)
Purity	>85% (SDS-PAGE)
Sequence	MCEYSKALHI LLKSPVTDDI IKFLTDTTLR VVPSSNYPTP PGSPGEKHLT RLPSLMTFIT RLVRYTNVYT PTLTAACYL NKLKRILPRD ATGLPSTIHR IFLACLILSA KFHNDSSPLN KHWARYTDGL FTLEDINLME RQLLQLLNWD LRVNTEDLIL DLQPLLEPIK QDLARSSDQR KRINMMMSMN RRTCAGTSP RSNNRFKLYE KQRNVSIA SD LSSATLV DSC NDLRRLK DVT NIANNTVANT NYVRTVEKWN DNVNRQSWDL EQIMSQHGF
Source	Mammalian cell
Target Names	PCL1
Protein Names	Recommended name: PHO85 cyclin-1 Alternative name(s): Cyclin HCS26 G1/S-specific cyclin PCL1
Expression Region	1-279
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.