




# Recombinant Drosophila melanogaster Ubiquitin-conjugating enzyme E2-17 kDa (UbcD6)

<b>Product Code</b>	CSB-YP328517DLU
<b>Relevance</b>	Catalyzes the covalent attachment of ubiquitin to other proteins. Required for postreplication repair of UV-damaged DNA.
<b>Abbreviation</b>	Recombinant Drosophila melanogaster UbcD6 protein
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	P25153
<b>Alias</b>	Ubiquitin carrier proteinUbiquitin-protein ligase
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Drosophila melanogaster (Fruit fly)
<b>Purity</b>	≥ 90% as determined by SDS-PAGE.
<b>Sequence</b>	MSTPARRRLMRDFKRLQEDPPTGVSGAPTENNIMIWNAVIFGPHDTPFEDGTF KLTIEFTEEYPNKPPTVRFVSKVFHPNVYADGGICLDILQNRWSPTYDVSAILTS IQLLSDPNPNSPANSTAAQLYKENRREYEKRVKACVEQSFID
<b>Research Area</b>	Others
<b>Source</b>	Yeast
<b>Target Names</b>	Ubc6
<b>Protein Names</b>	Recommended name: Ubiquitin-conjugating enzyme E2-17 kDa EC= 6.3.2.19 Alternative name(s): Ubiquitin carrier protein Ubiquitin-protein ligase
<b>Expression Region</b>	1-151aa
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	N-terminal 6xHis-tagged
<b>Mol. Weight</b>	19.2kDa
<b>Protein Length</b>	Full Length
<b>Image</b>	



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

### Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at  $-20^{\circ}\text{C}/-80^{\circ}\text{C}$ . Our default final concentration of glycerol is 50%. Customers could use it as reference.

### Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at  $-20^{\circ}\text{C}/-80^{\circ}\text{C}$ . The shelf life of lyophilized form is 12 months at  $-20^{\circ}\text{C}/-80^{\circ}\text{C}$ .