



Recombinant *Saccharomyces cerevisiae* Nucleosome assembly protein (NAP1)

Product Code	CSB-MP328528SVG
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P25293
Product Type	Recombinant Protein
Immunogen Species	<i>Saccharomyces cerevisiae</i> (strain ATCC 204508 / S288c) (Baker's yeast)
Purity	>85% (SDS-PAGE)
Sequence	MSDPIRTKPK SSMQIDNAPT PHNTPASVLN PSYLKNGNPV RAQAQEQQDDK IGTINEEDIL ANQPLLLQSI QDRLGSLVGQ DSGYVGGLPK NVKEKLLSLK TLQSELFVE KEFQVEMFEL ENKFLQKYKP IWEQRSRIIS GQE QPKPEQI AKGQEIVESL NETELLVDEE EKAQNDSEEE QVKGIPSFWL TALENLPIVC DTITDRDAEV LEYLQDIGLE YLTDGRPGFK LLFRFDSSAN PFFTNDILCK TYFYQKELGY SGDFIYDHA E GCEISWKDNA HNVTVDLEMR KQRNKTTKQV RTIEKITPIE SFFNFFDPPK IQNEDQDEEL EEDLEERLAL DYSIGEQ LKD KLIPRAVDWF TGAALFEFE EDEEEADEDE DEEEDDDHGL EDDDGESAE E QDDFAGRPEQ APECKQS
Source	Mammalian cell
Target Names	NAP1
Protein Names	Recommended name: Nucleosome assembly protein
Expression Region	1-417
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.