



Recombinant *Saccharomyces cerevisiae* DNA repair protein RAD57 (RAD57)

Product Code	CSB-MP328529SVG
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P25301
Product Type	Recombinant Protein
Immunogen Species	<i>Saccharomyces cerevisiae</i> (strain ATCC 204508 / S288c) (Baker's yeast)
Purity	>85% (SDS-PAGE)
Sequence	MPRALSIKFD NTYMDLYDEL PESKLLYDEE FSYLLDAVRQ NGVCVVDFLT LTPKELARLI QRSINEVFRF QQLLVHEYNE KYLEICEKNS ISPDNGPECF TTADVAMDEL LGGGIFTHGI TEIFGESSTG KSQLLMQLAL SVQLSEPAGG LGGKCVYITT EGDLPQRL E SMLSSRPAYE KLGITQSNIF TVSCNDLINQ EHIINVQLPI LLETSKGSIK LVIIDSISHH LRVELQNKSF RESQENKNYL DRMAEKLQIL AHDYSLSVVV ANQVGDKPLA NSPVAHRTYV TDYDYQLGWL VGWKNSTILY RQMNSLLGAS SNNDEILSDD EDYMLIERVM STVNDRNYDF FSKKKPPIIE NKTVERNSSS PISRQSKKRK FDYRVPNLGL TWSNHVSTRI LLQKSFKAST IQRGEAHL Y KGGDSASFWQ VKRTMKVVYS TFAKPGQIAY QITKRG IETA
Source	Mammalian cell
Target Names	RAD57
Protein Names	Recommended name: DNA repair protein RAD57
Expression Region	1-460
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.