



Recombinant Human Tubulin alpha chain-like 3 (TUBAL3)

Product Code	CSB-MP025315HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	A6NHL2
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MRECLSIHIG QAGIQIGDAC WELYCLEHGI QPNGVVLDTQ QDQLENAKME HTNASFDFFF CETRAGKHVP RALFVDLEPT VIDGIRTGQH RSLFHPEQLL SGKEDAANNY ARGRYSVGSE VIDLVLERTR KLAEQCGGLQ GFLIFRSFGG GTGSGFTSLL MERLTGEYSR KTKLEFSVYP APRISTAVVE PYNVLTTHS TTEHTDCTFM VDNEAVYDIC HRKLGVECPS HASINRLVVQ VVSSITASLR FEGPLNVDLI EFQTNLVPYP RIHFPMATAFA PIVSADKAYH EQFSVSDITT ACFESSNQLV KCDPRLGKYM ACCLLYRGDV VPKEVNAAIA ATKSRHSVQF VDWCPGTGFKV GINNRPPTVM PGGDLAKVHR SICMLSNTTA IVEAWARLDH KFDLMYAKRA FLHWYLREGM EEAFLARE DLAALERDYE EVAQSF
Source	Mammalian cell
Target Names	TUBAL3
Protein Names	Recommended name: Tubulin alpha chain-like 3
Expression Region	1-446
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.