



# Recombinant Escherichia coli Anaerobic sulfatase-maturing enzyme homolog AsIB (asIB)

<b>Product Code</b>	CSB-YP328549ENV
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	P25550
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Escherichia coli (strain K12)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MLQQVPTRAF HVMAKPSGSD CNLNCDYCFY LEKQSLYREK PVTHMDDDTL EAYVRHYIAA SEPQNEVAFT WQGGEPTLLG LAFYRRAVAL QAKYGAGRKI SNSFQTNGVL LDDEWCAFLA EHHFLVGLSL DGPPEIHNQY RVTKGGRPTH KLVMRALTLL QKHHVDYNVL VCVNRTSAQQ PLQVYDFLCD AGVEFIQFIP VVERLADETT ARDGLKLHAP GDIQGELTEW SVRPEEFGEF LVAIFDHWIK RDVGKIFVMN IEWAFANFVG APGAVCHHQP TCGRSVIVEH NGDVYACDHY VYPQYRLGNM HQQTIAEMID SPQQQAFGED KFKQLPAQCR SCNVLKACWG GCPKHRFMLD ASGKPLNYL CAGYQRYFRH LPPYLKAMAD LLAHGRPASD IMHAHLLVVS K
<b>Source</b>	Yeast
<b>Target Names</b>	asIB
<b>Protein Names</b>	Recommended name: Anaerobic sulfatase-maturing enzyme homolog AsIB Short name= AnSME homolog
<b>Expression Region</b>	1-411
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.