



Recombinant *Saccharomyces cerevisiae* Chromatin structure-remodeling complex protein RSC6 (RSC6)

Product Code	CSB-MP328554SVG
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P25632
Product Type	Recombinant Protein
Immunogen Species	<i>Saccharomyces cerevisiae</i> (strain ATCC 204508 / S288c) (Baker's yeast)
Purity	>85% (SDS-PAGE)
Sequence	MVTQTNPVPV TYPTDAYIPT YLPDDKVSNL ADLKKLIEMD SRLDLYLTRR RLDTSINLPT NTKTKDHPPN KEMLRIYVYN TTESSPRSDS GTPADSGKTT WTLRIEGKLL HESANGKHPF SEFLEGVAVD FKRLKPLGMG KKRKRDSLSL LPLNLQQPEY NDQDSTMGDN DNGEDEDSEAE AESREEIVDA LEWNYDENNV VEFDGIDIKR QGKDNLRCSI TIQLRGVDGG KVQYSPNLAT LIGMQTGSVN DAVYSIYKYI LINNLFVTEQ TEAQDGSNDA EDSSNENNNK NGAGDDDGVE GSTPKDKPEL GEVKLDSLLQ KVLDTNAAHL PLMNVVQTVN KLVSPLPPII LDYTIDLSKD TTYGATTLDV DVSHILHQPQ PQPNLQKEEE TDAEDTAKLR EITKLALQLN SSAQKYQFFH ELSLHPRETL THYLWSSKQN ELVLQGDQYF NEDAARTSDI YSNNNNDRSL MGNISLLYSQ GRL
Source	Mammalian cell
Target Names	RSC6
Protein Names	Recommended name: Chromatin structure-remodeling complex protein RSC6 Alternative name(s): Remodel the structure of chromatin complex subunit 6
Expression Region	1-483
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.