



# Recombinant *Saccharomyces cerevisiae* Mitotic check point protein BUB2 (BUB2)

<b>Product Code</b>	CSB-EP328620SVG
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	P26448
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	<i>Saccharomyces cerevisiae</i> (strain ATCC 204508 / S288c) (Baker's yeast)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MTSIEDLISN PLLLLHSSLS QLRYLILSEG LPISEDKQQQ RTRCYVWTVL SQTSMFASTQ RYLALLKLGP PSTTIYQKIK NDTSRFTQTD PNFRNRVSED ALIRCLSCFA WQTQRRRQKT RFGRIPVSTY VQGMNVLLAP LLYSCPSEPM AYQLFTKLCY EMIPTYLTKN LNGAQNGAKL LDISLRIIDP KLSKFLSDNL LTAEIYGMPIS ILTLSSCNKP LDQVIKLWDF MFAYGFHMNI LFVVAFLVKM RSKVFKSDSP VNLLRQFPDF DADEIIRLGV GFIAKIPAQI YDLLVDHLTD PDIYIP
<b>Source</b>	E.coli
<b>Target Names</b>	BUB2
<b>Protein Names</b>	Recommended name: Mitotic check point protein BUB2 Alternative name(s): Cell cycle arrest protein BUB2
<b>Expression Region</b>	1-306
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.