



Recombinant *Saccharomyces cerevisiae* 60S ribosomal protein L16-B (RPL16B)

Product Code	CSB-EP328651SVG
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P26785
Product Type	Recombinant Protein
Immunogen Species	<i>Saccharomyces cerevisiae</i> (strain ATCC 204508 / S288c) (Baker's yeast)
Purity	>85% (SDS-PAGE)
Sequence	SQPVVVIDA KDHLLGRLAS TIAKQVLNGQ KIVVVRAEAL NISGEFFRNK LK YHDFLRKA TAFNKTRGPF HFRAPSRILY KAIRGMVSHK TARGKAALER LKIFEGIPPP YDKKKRVVVP QALRVLR LKP GRKYTTLGKL STSVGWKYED VVAKLEDK RK VRS AEYYAKK RAFTKKVSSA SAAASESDVA KQLASFGY
Source	E.coli
Target Names	RPL16B
Protein Names	Recommended name: 60S ribosomal protein L16-B Alternative name(s): L21 RP23 YL15
Expression Region	2-198
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.