



Recombinant Type-2 restriction enzyme A_{br}I (abrIR)

Product Code	CSB-EP328662DPS
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P26919
Product Type	Recombinant Protein
Immunogen Species	Azospirillum brasilense
Purity	≥85% (SDS-PAGE)
Sequence	MALELADYER KARESVMAFW GNREKARQKQ QESGKPDQGE RRRDRRKNMD GFCALVIDII RANGLHRAEI HQKRALLTLP GYFRPTKLWD LLVIQDRRLV AALEFKSQVG PSFGNMFNNR TEEAIGTAHD LWTAYREGAF GQSPRPFVGW LMLVEDAPES NQRRADTSPH FRCSSRAH PTYGGMTSSA SGSSRSSFTP RPRSSPRRAP PPRQGRIISP RI
Source	E.coli
Target Names	abrIR
Protein Names	Recommended name: Type-2 restriction enzyme A _{br} I Short name= R.A _{br} I EC= 3.1.21.4 Alternative name(s): Endonuclease A _{br} I Type II restriction enzyme A _{br} I
Expression Region	1-222
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.