



Recombinant Crimean-Congo hemorrhagic fever virus Nucleoprotein (N)

Product Code	CSB-YP328701CSB
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P27317
Product Type	Recombinant Protein
Immunogen Species	Crimean-Congo hemorrhagic fever virus (isolate C68031) (CCHFV)
Purity	>85% (SDS-PAGE)
Sequence	MENKIEVNNK DEMNKNWFEEF KKGNGLVDTF TNPYSFCESV PNLERFVFQM ASATDDAQKD SIYASALVEA TKFCAPIYEC AWWSSSTGIVK KGLEWFEKNA GTIKSWDESY IELKVEVPKI EQLANYQQAA LKWRKDIGFR VNANTAALSH KVLAEYKVPK EIVMSVKEML SDMIRRRNLI LNRGGDENPR GPVSREHVEW CREFVKGKYE MAFNPPWGGI NKSGRSGIAL VATGLAKLAE TEGKGVFDEA KKTVEALNGY LDKHKDEVDK ASADNMITNL LKHIAKAQEL YKNSSALRAQ GAQIDTAFSS YYWLYKAGVT PETFPTVSQF LFELGKQPRG TKMKMKALLS TPMKWGKKLY ELFADDSFQQ NRIYMHPAVL TAGRISEMGV CFGTIPVANP DDAAQGSQHT KSILNLRNT ETNNPCAkti VKLFEIQKTG FNIQDMDIVA SEHLLHQSLV GKQSPFQNAV NVKGNATSAN II
Source	Yeast
Target Names	N
Protein Names	Recommended name: Nucleoprotein Alternative name(s): Nucleocapsid protein Short name= Protein N
Expression Region	1-482
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.