



Recombinant Renilla reniformis Coelenterazine h 2-monooxygenase

Product Code	CSB-MP328734RCE
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P27652
Product Type	Recombinant Protein
Immunogen Species	Renilla reniformis (Sea pansy)
Purity	>85% (SDS-PAGE)
Sequence	MTSKVYDPEQ RKRMITGPQW WARCKQMNVL DSFINYYDSE KHAENAVIFL HGNAASSYLW RHHVPHIEPV ARCIIPDLIG MGKSGKSGNG SYRLLDHYKY LTAWFELLNL PKKIIFVGHD WGACLAFHYS YEHQDKIKAI VHAESVVDVI ESWDEWPDIE EDIALIKSEE GEKMLENNF FVETMLPSKI MRKLEPEEFA AYLEPFKEKG EVRRPTLSWP REIPLVKGGK PDVVQIVRNY NAYLRASDDL PKMFIESDPG FFSNAIVEGA KKFPNTEFVK VKGLHFSQED APDEMGKYIK SFVERVLKNE Q
Source	Mammalian cell
Protein Names	Recommended name: Renilla-luciferin 2-monooxygenase EC= 1.13.12.5Alternative name(s): Renilla-type luciferase
Expression Region	1-311
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.