



# Recombinant *Oryza sativa* subsp. *japonica* Alpha-amylase isozyme 3D (AMY1.3)

<b>Product Code</b>	CSB-EP328760OFG
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	P27933
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	<i>Oryza sativa</i> subsp. <i>japonica</i> (Rice)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	QVLFQ GFNWESWKQQ GGWYNMLKGQ VDDIAKAGVT HVWLPPPSHS VAPQGYMPGR LYDL DASKYG TAAELKSLIA AFHGKGVQCV ADVVINHRCA EKKDARGVYC VFEGGTPDDR LDWGPGMICS DDTQYSDGTG HRDTGEGFGA APDIDHLNPR VQRELTDWLN WLKSDVGF DG WRLDFAKGYS TDIKMYVES CKPGFVVAEI WNSLSYNGDG KPAANQDQGR QELVNWVNAV GGPAMTFDFT TKGLLQAGVQ GELWRLRDGN GKAAGMIGWL PEKAVTFVDN HDTGSTQKLW PFP SDKVMQG YAYILTHPGV PCIFYDHMFD WNLKQEITAL AAIRERNGIN AGSKLRIVVA DADAYVAVVD EKVMVKIGTR YDVGNAVPSD FHQTVHGKDY SVWEKGS LRV PAGRHL
<b>Source</b>	E.coli
<b>Target Names</b>	AMY1.3
<b>Protein Names</b>	Recommended name: Alpha-amylase isozyme 3D EC= 3.2.1.1 Alternative name(s): 1,4-alpha-D-glucan glucanohydrolase
<b>Expression Region</b>	26-436
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full Length of Mature Protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.