



Recombinant *Saccharomyces cerevisiae* SWIRM domain-containing protein FUN19 (FUN19)

Product Code	CSB-YP328767SVG
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P28003
Product Type	Recombinant Protein
Immunogen Species	<i>Saccharomyces cerevisiae</i> (strain ATCC 204508 / S288c) (Baker's yeast)
Purity	>85% (SDS-PAGE)
Sequence	MGLYSPSESEK SQLNMNYIGK DDSQSIFRRL NQNLKASNNN NDSNKNGLNM SDYSNNSPYG RSYDVRINQN SQNNGNGCFS GSIDSLVDEH IIPSPPLSPK LESKISHNGS PRMASSVLVG STPKGAVENV LFVKPVWPNG LSRKRYRYAT YGFLSQYKIF SNLAQPYSKN IINRYNNLAY NARHKYSKYND DMTTPPLPS SSSRLPSPLA SPNLNRQARY NMRKQALYNN NLGKFESDTE WIPRKRKVYS PQRRTMTTSP HRAKKFSPSA STPHTNIASI EAIHDAPQYI PNVSWKKLPD YSPPLSTLPT DSNKSLKIEW KGSPMDLSTD PLRNELHPAE LVLAQTLRLP CDLYLDSKRR LFLEKVVYRLK KGLPFRRRTDA QKACRIDV NK ASRLFQAF EK VGWLQDSNFT KYL
Source	Yeast
Target Names	FUN19
Protein Names	Recommended name: SWIRM domain-containing protein FUN19
Expression Region	1-413
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.