



Recombinant Human Nuclear pore complex-interacting protein family member B15 (NPIP-B15)

Product Code	CSB-EP015992HU-B
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	A6NHN6
Storage Buffer	Lyophilized from Tris/PBS-based buffer, 6% Trehalose, pH 8.0
Product Type	Recombinant Proteins
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	QP TPVINS LAVY RHRE TDFGVG VRDHPGQH GK TPSPQKLDNL IIIIGFLRR DTFTILFCTS YLCVSFLKTI FWSRNGHDGS TDVQQRAWRS NRSRQKGLRS ICMHTKKRVS SFRGNKIGLK DVITLRRHVE TKVRAKIRKR KVTTKINRHD KINGKRKTAR KQKMFQRAQE LRRRAEDYHK CKIPPSARKP LCNWVRMAAA EHRHSSGLPC WPYLTAEALK NRMGRQPPPP TQQHSITDNS LSLKTPPECL LHPLPPSVDD NIKCPLAPL PPSVDDNLKE YLLVPLPPSP LPPSVDDNLK DCLFVPLPPS PLPPSVDDNL KTPPLATQEA EAEKPPKPKR WRVDEVEQSP KPKRRRADEV EQSPKPKRQR EAEAQQLPKP KRRRLSKLRT RHCTQAWAIR INPWVEKKKK IKKQNKTHAP KTN
Source	E.coli
Target Names	NPIP-B15
Protein Names	Recommended name: Nuclear pore complex-interacting protein-like 2
Expression Region	19-443
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.