



Recombinant *Candida albicans* Zinc finger protein 1 (CZF1)

Product Code	CSB-MP328848CZE
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P28875
Product Type	Recombinant Protein
Immunogen Species	<i>Candida albicans</i> (strain WO-1) (Yeast)
Purity	≥85% (SDS-PAGE)
Sequence	MSSIPNINWN DPNNGKSNTS RQSQPQPQLP SNVSPSNSRA VPTSGSIGGP QYGSSQFSNE YSRNPNTIGG PPFPLQSNQR GYPNTGYPV QQTAQQRSGD KLQQVHSQQQ QQQQQQQQPL YQQYPPQSVG YLAGDVYNPQ HQEYVQMNQL PNQHYNLQQR QQAQGGQLKS QLNEQNAMMS ASTQQYPVQD FTNPYPNAQN PAEQQQQQQP LRTQSQQWDG YQSQPLYSAA GNTIPSSIQQ QIPPQNLSPS EQQQVKQQQP SPPEQGTKKK PGRKPKLRKL SESSSETPQV PKTASSSSSS PTAVNSGKPI TKRSRMGCLT CRQRKKRCCE TRPRCTECTR LRLNCTWPKP GTEHKNPKD QKDDENTIEH AEFGRIVLR GIVEYRSK
Source	Mammalian cell
Target Names	CZF1
Protein Names	Recommended name: Zinc finger protein 1
Expression Region	1-388
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.