



Recombinant Acyl-CoA reductase (luxC)

Product Code	CSB-YP328885EUW
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P29236
Product Type	Recombinant Protein
Immunogen Species	Photobacterium leiognathi
Purity	≥85% (SDS-PAGE)
Sequence	MIKKIPMIIG GVVQNTSGYG MRELTLNNNK VNIPITQSD VEAIQSLNIE NKLTIQIVN FLYTVGQKWK SETYSRRLTY IRDLIKFLGY SQEMAKLEAN WISMILCSKS ALYDIVENDL SSRHIIDEWI PQGECYVKAL PKGKSVHLLA GNVPLSGVTS ILRAILTKNE CIIKTSSADP FTATALVNSF IDVDAEHPIT RSISVMYWSH SEDLAIPKQI MSCADVVIW GGDDAIKWAT EHAPSHADIL KFGPKKSISI VDNPTDIKAA AIGVAHDICF YDQQACFSTQ DIYYIGDSID IFFDELAQQL NKYKDILPKG ERNFDEKAAF SLTERECLFA KYKVQKGESQ SWLLTQSPAG SFGNQPLSRS AYIHQVNDIS EVIPFVHKAV TQTVAIAPWE SSFKYRDILA EHGAERIEA GMNNIFRVGG AHDGMRPLQR LVNYISHERP STYTTKDVSV KIEQTRYLEE DKFLVFVP
Source	Yeast
Target Names	luxC
Protein Names	Recommended name: Acyl-CoA reductase EC= 1.2.1.50
Expression Region	1-478
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.