



Recombinant Xylose isomerase (xylA)

Product Code	CSB-EP328904KBG
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P29442
Product Type	Recombinant Protein
Immunogen Species	Klebsiella pneumoniae
Purity	>85% (SDS-PAGE)
Sequence	<pre>MQTYFDQLDR VRYEGPKSTN PLAFRHYNPD ELVLGKRMED HLRFAACYWH TFCWNGADMV GVGSFDRPWQ QPGDALEMAK RKADVAFEFF HKLNVPIYCF HDVDVSPEGA SLKEYSNNFA RMVEVLAEKQ QQSGVKLLWG TANCFTNPRY GAGAATNPDP EVFSWAATQV VTAMNATHQL GGENYVLWGG REGYETLLNT DLRQEREQLG RLLQLVVEHK HKIGFKGTLL IEPKPQEPTK HQYDYDASTV YGFLKQFGL EKEIKLNIEAN HATLAGHSFH HEIATAIALG LFGSVDANRG DPQLGWDTDQ FPNRVEEDAL VMYEIFKAGG FTTGGLNFDA KVRREQSTDKY DLFYGHIGAM DTMAVSLKVA ARMIEDGELD KRVARRYAGW NGELGQQILN GQMTLSDIAQ YAAQHQLAPQ HRSGQQEQLE NLVNHLYLFDK</pre>
Source	E.coli
Target Names	xylA
Protein Names	Recommended name: Xylose isomerase Short name= XI EC= 5.3.1.5
Expression Region	1-440
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.