



# Recombinant Vaccinia virus Major envelope protein (F5)

<b>Product Code</b>	CSB-MP328938VAD
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	P29885
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Vaccinia virus (strain L-IVP) (VACV)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MWPFASVPAG AKCRLVETLP ENMDFRSDHL TTFECFNEII TLAKKYIYIA SFCCNPLSTT RGALIFDKLK EASEKGIIK VLLDERGKRN LGELQSHCPD INFITVNIDK KNNVGLLLGC FWVSDDERCY VGNASFTGGS IHTIKTLGVY SDYPPLATDL RRRFDTFKAF DSAKNSWLNL CSAACCLPVS TAYHIKNPIG GVFFTDSEPH LLGYSRDLDT DVVIDKLSA KTSIDIEHLA IVPTRVDGN SYYWPDIYNS IIEAAINRGV KIRLLVGNWD KNDVYSMATA RSLDALCVQN DLSVKVFTIQ NNTKLLIVDD EYVHITSANF DGTHYQNHGF VSFNSIDKQL VSEAKKIFER DWVSSHSLKI
<b>Source</b>	Mammalian cell
<b>Target Names</b>	F5
<b>Protein Names</b>	Recommended name: Major envelope protein Alternative name(s): 37 kDa protein Palmitoylated EEV membrane protein p37K
<b>Expression Region</b>	1-372
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.