



Recombinant Solanum tuberosum Threonine dehydratase biosynthetic (TD)

Product Code	CSB-MP329054FIG
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P31212
Product Type	Recombinant Protein
Immunogen Species	Solanum tuberosum (Potato)
Purity	>85% (SDS-PAGE)
Sequence	PFDAPGVKIG QGTIGTEINR QLKDIHAVFV PVGGGGLISG VAAYFTQVAP HTKIIGVEPY GAASMTLSLY EGHRVKLENV DTFADGVAVA LVGEYTFACK QELIDGMVLV RNDGISAAIK DVYDEGRNIL ETSGAVAIAG AAAYCEFYN KNENIVAIAS GANMDFSCLKH KVTELAELGS DNEALLATFM IEQPGSFKTF AKLVGSMNIT EVTYRFTSER KEALVLYRVD VDEKSDLEEM IKKLNSSNMK TFNFHSHNELV AEHIKHLVGG SASISDEIFG EFIFPEKAGT LSTFLEAFSP RWNITLCRYR DQGDINGNVL VGFQVPQSEM DEFKSQADGL GYPYELDNSN EAFNIVVAE
Source	Mammalian cell
Target Names	TD
Protein Names	Recommended name: Threonine dehydratase biosynthetic EC= 4.3.1.19 Alternative name(s): Threonine deaminase Short name= TD
Expression Region	1-359
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.